

Pyruvate Assay Kit

Catalog # EA-7020

Introduction

The Pyruvate Assay Kit utilizes an enzyme catalyzed oxidation reaction to measure pyruvate levels in samples. First, pyruvate is oxidized by pyruvate oxidase to form a pyruvate byproduct and hydrogen peroxide. Then, the pyruvate level in the sample is determined by quantifying the hydrogen peroxide generated by the enzyme reaction with a fluorogenic probe that can be measured with a spectrophotometer.

Materials Required but Not Provided

- PBS
- 96-well clear microplate for absorbance reading or 96well black microplate with clear bottom for fluorescence reading
- Microplate reader capable of measuring absorbance at 560 nm or fluorescence at 530nm/590nm

Materials Provided

- 10mM FAD (-80°C)
- Probe Reagent (-20°C)
- HRP Reagent (4°C)
- 20mM Pyruvate Standard (-20°C)
- POX Enzyme Stock (-80°C)

**Spin down small tubes before starting experiment. **

(For Research Use Only)

Plasma Sample Preparation

- 1. Centrifuge citrated or EDTA-collected blood at 4° C (1,000 x g for 10 minutes) to separate plasma from erythrocytes. Alternatively, blood collected without anticoagulant can be centrifuged to collect serum.
- 2. Transfer the plasma layer to a new tube without disturbing the buffy layer.
- 3. The plasma may be assayed directly or stored away at -80°C.

Cell Sample Preparation

- 1. Wash the cells once with PBS before lysing the cells.
- 2. For a 96-well culture plate, add 40 μL of Lysis buffer to each well and incubate at room temperature for 10 minutes.
- 3. Pipette the Lysis buffer up and down to detach the cells and transfer the cell lysates into a new tube.
- 4. If necessary, homogenize the cell lysates with a sonicator.
- 5. The cell lysates may be assayed directly or stored at -80° C.

Tissue Sample Preparation

- 1. Weigh tissue sample and add 1 mL of Tissue Lysis buffer per 100mg of tissue.
- 2. Homogenize the tissue samples with a tissue grinder.
- 3. If necessary, further homogenize the tissue samples with a sonicator.
- 4. Centrifuge the sample at 10,000 RPM for 5 minutes to pellet the tissue debris.
- 5. Collect the supernatant and measure the protein concentration of the supernatant. The tissue sample can be assayed directly or stored at -80°C.
- 6. Use the Dilution buffer to dilute the tissue sample to the appropriate concentration for each assay.

support@signosisinc.com Technical Support 1(408)-747-0771 Telephone

Pyruvate Measurement

1. <u>Standard curve preparation</u>: First, take eight new tubes and label them 1-8. In tube 1, prepare a 2000 μ M pyruvate standard by mixing 10 μ L of the provided 20 mM pyruvate standard with 90 μ L of PBS. Then, for tubes 2-8, add 50 μ L of PBS to each tube. Transfer 50 μ L of the 2000 μ M pyruvate standard from tube 1 to the PBS in tube 2 and mix to make a 1000 μ M pyruvate standard. Then, transfer 50 μ L of the 1000 μ M pyruvate standard from tube 2 to the PBS in tube 3 and mix to make a 500 μ M pyruvate standard. Continue the serial dilutions until tube 7 is done. Leave tube 8 untouched as the PBS negative control.

Standard#	Pyruvate
	Concentration (µM)
1	2000
2	1000
3	500
4	250
5	125
6	62.5
7	31.25
8	0

2. <u>Reaction mix preparation</u>: calculate the amount of each reagent needed to make the reaction mix according to the table below.

Component	Reaction Mix (per well/sample)
10mM FAD	0.05 μL
POX Enzyme	0.05 μL
PBS	49.9 μL
Total	50 µL

- 3. Any unused enzyme stock can be stored at -80°C for future use.
- 4. Add 50 μ L of reaction mix to each well of the plate.
- 5. Add 5 μ L of sample or standard to each well with reaction mix and mix thoroughly.
- 6. Cover the plate and incubate at room temperature for 30 minutes.
- 7. <u>Detection mix preparation</u>: calculate the amount of each reagent needed to make the detection mix according to the table below.

Component	Detection Mix (per well/sample)
Probe Reagent	0.5 μL
HRP Reagent	1 μL
PBS	48.5 μL
Total	50 µL

www.signosisinc.com Home Page info@signosisinc.com Questions / Comments

- 8. Add 50 μ L of detection mix to each reaction well in the plate. Be sure to add the detection mix quickly, since the signal begins to develop when the reagents are added. Use a multichannel pipette if possible.
- Cover the plate and incubate at room temperature away from light for 10-20 minutes.
 Exposure to light will produce background signal in wells
- 10. For a stronger signal, the plate can be incubated for another 30-60 minutes away from light.
- 11. Measure the absorbance of the plate at 560 nm using a plate reader. Alternatively, measure the fluorescence of the plate in a fluorescence plate reader Ex/Em 530nm/590nm.

support@signosisinc.com Technical Support 1(408)-747-0771 Telephone